

مؤسسة المهند للأجهزة الطبية  
AL MOHANAD MEDICAL EQUIPMENT EST.



جودة أفضل لحياة أفضل  
Good Quality for Better Life

تقرير الصيانة  
Maintenance Report

No. 02698

Work Order No.: 1  
Equipment: Accent 200  
Manufacturer: Corning Diagnostic  
Machine No.: 1

Date: 05/12/2022  
Location: Olas med. center  
Model: Chem  
Serial No.: WVN-4C105822

Eng's. Action

Requested By:

checked the device. Replaced the filter. Found issue in cooling module. The results stability is due to cooling failure not good. Also found creation of B/T in refrigerator should be in room temp. Please follow the standard op to increase the stability of reagents. Tared the staff on how to use device & load values.

The following Parts are

☒ Required

☐ Changed

Part #	Description	Qty.
-	Pelton	2
	Recommend: waste tube, DI tube	

Job Completed on 05/12/2022 At am/pm

Repair Time: From to

Travel Time: From to

Eng. Name / Sig.

☐ Work Completed, Unit Left Operational.

☒ Work Un Completed, Parts Required.

This Jop is Done

Clinic / Lab Incharge

Name & Date:

Approved By Eng.

